

ARTICLE INFORMATION SHEET

Page 1/2

Issuing Date 03/11/2021

Revision Date 03/11/2021

Revision Number 1.14

1. IDENTIFICATION OF THE ARTICLES

Product name DRYVIEW DVE Laser Imaging Film

Product Code(s) 4674

Recommended Use Medical X-ray Film.

2. COMPANY

Supplier: Carestream Health, Inc., 150 Verona Street, Rochester, NY, USA 14608

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical name	CAS No.	Weight-%
Polyethylene terephthalate film base	-	90-95
Polyvinyl butyral-polyvinyl alcohol-polyvinyl acetate terpolymer	-	1-5
Docosanoic acid, silver(1+) salt	-	1-5
Cellulose, acetate butanoate	-	<3

4. HEALTH AND SAFETY	
Potential Health Effects	PET films are not hazardous to health. No personal protection is needed under recommended conditions of use. May be harmful if swallowed.
Fire Hazard	As the product contains combustible organic components, fire will produce dense black smoke containing hazardous products of combustion. In the event of a fire, fire-fighters should wear self contained breathing apparatus and protective clothing.
Suitable Extinguishing Media	Water spray. Carbon dioxide (CO2). Dry chemical.
5. WASTE DISPOSAL	

Dispose of in accordance with local regulations. Recover silver before disposal.

6. TRANSPORT AND LABELING REGULATIONS

DOT	Not regulated
TDG	Not regulated
IATA	Not regulated
IMDG	Not regulated
ADR	Not regulated

7. STORAGE

Keep in a dry, cool place.

8. OTHER INFORMATION

The information contained herein is furnished without warranty of any kind. Users should consider these data only as a supplement to other information gathered by them and must make independent determinations of suitability and completeness of information from all sources to assure proper use and disposal of these materials and the safety and health of employees and customers and the protection of the environment.