

ARTICLE INFORMATION SHEET

Page 1 / 2

Issuing Date 2015-01-08

Revision date 2015-01-08

Version 1.01

1. IDENTIFICATION OF THE ARTICLE

Product Name: DRYVIEW DVB+

Product Code: 4678

Identified uses: Photographic film. For industrial use only.

2. COMPANY

Supplier Carestream Health, Inc., 150 Verona Street, Rochester, New York 14608

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS-No	Weight percent
Polyethylene terephthalate film base	-	90-95
Polyvinyl butyral-polyvinyl alcohol-polyvinyl acetate terpolymer	-	1-5
Docosanoic acid, silver(1+) salt	-	1-5
Cellulose, acetate butanoate	-	<3

4. HEALTH AND SAFETY

Suitable extinguishing media Water spray, Carbon dioxide (CO₂), Dry chemical

5. WASTE DISPOSAL

Dispose of in accordance with local regulations.

6. TRANSPORT AND LABELLING REGULATIONS

DOT Not regulated

TDG Not regulated

ICAO/IATA Not regulated

IMDG/IMO Not regulated

ADR/RID Not regulated

For transportation information, go to: <http://ship.carestream.com>

7. STORAGE

Keep in a dry, cool place.

8. OTHER INFORMATION

The information contained herein is furnished without warranty of any kind. Users should consider the data only as a supplement to other information gathered by them and must determine independently suitability and completeness of information from all sources to ensure proper use and disposal of these materials, the safety and health of employees and customers and the protection of the environment.